

**DTU Food**  
National Food Institute



**PROTOCOL FOR WHOLE GENOME SEQUENCING AND BIOINFORMATIC ANALYSIS OF BACTERIAL ISOLATES RELATED TO THE EU MONITORING OF ANTIMICROBIAL RESISTANCE**

**AUTHORED BY THE EURL-AR**

Jette Sejer Kjeldgaard

[jetk@food.dtu.dk](mailto:jetk@food.dtu.dk)

National Food Institute, DTU  
European Union Reference Lab on  
Antimicrobial Resistance

<https://www.eurl-ar.eu/wgs.aspx>

# Background for protocol

- The Commission Implementing **Decision 2020/1729** on the monitoring and reporting of antimicrobial resistance in zoonotic and commensal bacteria
- Authorising the **use of WGS as an alternative method** for prediction of resistance in relation to the specific monitoring of ESBL- or AmpC- or carbapenemase-producing *E. coli* and *Salmonella*
- The EURL-AR has produced the present protocol for guidance in these matters
- The whole genome sequencing (WGS) processes divides into three overall processes:
  - Bacterial isolation, DNA preparation and DNA quality and quantity assessment
  - Library preparation, library quality and quantity assessment and sequencing
  - Sequence QC and bioinformatics analyses

# Purpose of protocol

- **Ensure that WGS data reported to EFSA is obtained in a harmonised and comparable way**
  - **Less important**
    - How the bacteria, DNA and sequences are obtained
  - **Very important**
    - Assure the sequence quality control
      - Using the same QC criteria
    - Harmonised AMR gene analysis
      - using the same methods and settings for analysis
    - Reporting adequate data

# Table of Contents

Important notes .....	4
Protocol .....	5
Bacterial isolation .....	5
DNA preparation and quality assessment .....	5
Library preparation .....	6
Sequencing .....	6
Assessment of the genomic sequence quality .....	7
AMR gene and point mutation prediction .....	9
Additional analysis and sub-typing .....	10
Proficiency test .....	10
Online training .....	10
References .....	11
Abbreviations and acronyms .....	12
Links .....	13

# Links Generic protocol – not one method that fits all

Table 2: Collection of links referred to in the protocol, including last date of accession

Link#	Method or content	Last accessed
Link 1	Illumina website <a href="https://www.illumina.com/">https://www.illumina.com/</a>	December 2020
Link 2	Oxford Nanopore website <a href="https://nanoporetech.com/products">https://nanoporetech.com/products</a>	December 2020
Link 3	Thermofisher website <a href="https://www.thermofisher.com/dk/en/home/life-science/sequencing/next-generation-sequencing/ion-torrent-next-generation-sequencing-products-services.html">https://www.thermofisher.com/dk/en/home/life-science/sequencing/next-generation-sequencing/ion-torrent-next-generation-sequencing-products-services.html</a>	December 2020
Link 4	EURL-AR website – Inter-EURLs WG on NGS <a href="https://www.eurl-ar.eu/inter-eurls-working-group-on-ngs.aspx">https://www.eurl-ar.eu/inter-eurls-working-group-on-ngs.aspx</a>	December 2020
Link 5	Document on bioinformatics tools for basic analysis of Next Generation Sequencing data <a href="https://www.iss.it/documents/20126/0/Bioinformatics_tools_for_basic_analysis_of_Next_Generation_Sequencing_data_Del4.pdf/02c8f77b-db2c-6b8d-e2ba-416144f89f7e?t=1602603602556">https://www.iss.it/documents/20126/0/Bioinformatics_tools_for_basic_analysis_of_Next_Generation_Sequencing_data_Del4.pdf/02c8f77b-db2c-6b8d-e2ba-416144f89f7e?t=1602603602556</a>	December 2020
Link 6	Methods for isolation of ESBL, ampC and carbapenemase-producing E. coli from meat and caecal samples <a href="https://www.eurl-ar.eu/protocols.aspx">https://www.eurl-ar.eu/protocols.aspx</a>	December 2020
Link 7	Method for detection of <i>Salmonella</i> in food and animal feed <a href="https://www.eurlsalmonella.eu/publications/eurl-manual">https://www.eurlsalmonella.eu/publications/eurl-manual</a>	December 2020
Link 8	Method for detection of Campylobacter <a href="https://www.sva.se/en/about-us/eurl-campylobacter/laboratory-procedures/">https://www.sva.se/en/about-us/eurl-campylobacter/laboratory-procedures/</a>	December 2020
Link 9	DNA extraction protocol EasyDNA <a href="https://assets.thermofisher.com/TFS-Assets/LSG/manuals/easydna_man.pdf">https://assets.thermofisher.com/TFS-Assets/LSG/manuals/easydna_man.pdf</a>	December 2020
Link 10	Automated DNA extraction Magna Pure <a href="https://lifescience.roche.com/en_dk/products/magna-pure-96-instrument-382411-1.html">https://lifescience.roche.com/en_dk/products/magna-pure-96-instrument-382411-1.html</a>	December 2020
Link 11	Overview of applications of Qubit <a href="http://www.invitrogen.com/qubit">www.invitrogen.com/qubit</a>	December 2020
Link 12	Protocol for Qubit 4 DNA quantification <a href="https://assets.thermofisher.com/TFS-Assets/BID/manuals/MAN0017210_Qubit_4_Assays_QR.pdf">https://assets.thermofisher.com/TFS-Assets/BID/manuals/MAN0017210_Qubit_4_Assays_QR.pdf</a>	December 2020

## Links to protocols

- Continuously updated
- Dependent on lab
  - Equipment
  - Throughput
  - Prerequisites

# Table of Contents

Important notes .....	4
Protocol .....	5
Bacterial isolation .....	5
DNA preparation and quality assessment .....	5
Library preparation .....	6
Sequencing .....	6
Assessment of the genomic sequence quality .....	7
AMR gene and point mutation prediction .....	9
Additional analysis and sub-typing .....	10
Proficiency test .....	10
Online training .....	10
References .....	11
Abbreviations and acronyms .....	12
Links .....	13

# Bacterial isolation, DNA preparation and DNA quality and quantity assessment

- Methods for isolation of ESBL, ampC and carbapenemase-producing E. coli from meat and caecal samples (EURL-AR)
- Method for detection of Salmonella in food and animal feed (EURL-Salmonella)
  - Ensure purity and correct species
- Examples of DNA extraction kits
  - Laboratory routine methods
- Examples of DNA quality/quantity assessment

# Library preparation, library quality and quantity assessment and sequencing

- Dependent on the laboratory equipment
  - Majority using Illumina sequencing equipment

At present the EURL-AR recommends Illumina sequencing

- QC of sequences
- Tools for analysis

- Suggestion for library preparation
- Quantification and QC of library prep
- Illumina instrument-specific sequencing reagents, flow cells, cluster generation reagents
  - MiSeq and NextSeq



# Sequence QC and bioinformatics analyses

- Trimming of raw reads
  - Can be performed, but is not crucial for Illumina sequences
- File format
  - it is recommended to perform the assembly of fastq files into fasta files
    - part of the quality control
- Check for contamination
  - E.g. using KmerFinder for species determination and look into QC parameters
- Assembly
  - Using SPAdes 3.14 or newer
  - Accessible as CGE tool with output of important QC parameters
    - <https://cge.cbs.dtu.dk/services/SPAdes-3.14/>

# QC parameters

- The process of raw reads assembly into contigs outputs a range of QC parameters
- **number of reads**
- **depth of coverage**
- **average read length** (as specified by the sequencing equipment)
- **size of assembled genome** (+/- 0.5 million bases deviation from expected size)
- **total number of contigs** (<500 contigs)
- **N50** (>30.000 bp)

# Assembly with SPAdes v 3.14

- The SPAdes 3.14 tool will output the contigs file (.fasta) and additionally a .txt file with some basic statistics and QC parameters.
- The output file contains data on:
  - Input files :
    - Total number of reads
    - Total number of bases
  - Contigs file :
    - Number of contigs
    - Number of bases (assembled genome size)
    - N50
- Using this output, it is also possible to calculate the average read length=  $\frac{\text{Number of bases}}{\text{Number of reads}}$  (input files)

# AMR gene and point mutation prediction

- The EURL-AR recommends using ResFinder v4.1 or newer
- For harmonisation of the AMR data reported by different laboratories, it is important to use the defined settings.
- The EURL-AR recommends running the ResFinder analysis on the contigs **assembly files (.fasta)** using specific **settings**
- ResFinder can be run as a web-tool (CGE) or as local installation (available on BitBucket)
  - Web-tool limited to analysing one sequence at a time

# ResFinder settings

For chromosomal point mutations:

- Select threshold for % ID: **90 %**
- Select minimum length: **60 %**

For acquired antimicrobial resistance genes:

Select all antimicrobial databases (default setting)

- Select threshold for % ID: **90 %**
- Select minimum length: **60 %**

Select species: as appropriate

Select type of your reads: Assembled genome/Contigs

**Chromosomal point mutations**

Select threshold for %ID  
90 %

Select minimum length  
60 %

Show unknown mutations, not found in the database

---

**Acquired antimicrobial resistance genes**

**Select Antimicrobial configuration**  
Select multiple items, with Ctrl-Click (or Cmd-Click on Mac) - as default a

Aminoglycoside  
Beta-lactam  
Colistin  
Fluoroquinolone  
Fosfomycin  
Fusidic Acid

Select threshold for %ID  
90 %

Select minimum length  
60 %

Acquired disinfectant resistance genes

---

**Select species**  
Campylobacter spp.\*  
\*Chromosomal point mutation database exists

**Select type of your reads**  
Assembled Genome/Contigs

# Data to report to EFSA

- Beyond the sampling and isolate data, the results reported in relation to Decision 2020/1729 should include:
  - Date of sequencing
  - Sequencing technology used
  - Library preparation used
  - Version of the predictive tool (ResFinder)
  - AMR-conferring genes data:
    - Gene name
    - Output information on % identity
    - Output information on % coverage (length)
  - Date of ResFinder analysis
- The protocol will be added a template sheet for collection of metadata, including examples of how to report data.

# Questions and discussion?